

# Detecting Brain Hemorrhage

Ramin Ostad Flatiron Data Science Fellowship




## **Problem**

"Hemorrhage or no hemorrhage?" That is the question...








## **Problem**

"Hemorrhage or no hemorrhage?" That is the question...







## Approach

We will train a convolutional neural network to identify signs of brain hemorrhage using 600K+ previously labeled brain scanned images

Image Preprocessing Convolutional
Neural Network

Prediction

## **Data**

- CT Scan images enter the model as DICOM files
- We extract a density level that allows our model to 'see' brain, blood and bones

## Model

- Model will be built in tensorflow using the Keras wrapper
- We will use various layers to detect similarities across same-class images

## **Hemorrhage Detection**

 The model can be used to predict whether a new patient has brain hemorrhage

## CT Scans of the Brain

CT scans combine many individual x-rays to produce a more detailed image of a specific area of the body

#### What is a CT scan?

A Computed Tomography scan, otherwise known as a CT or CAT scan, is a process whereby many x-rays are taken of an body part from different angles and combined to form a single more detailed image of that body part

### How are CT scans taken?

CT scans are taken using a CT scanner. Typically a CT scanner consists of a tube which houses the x-ray detector. The x-ray detector moves around the body in the tube as the body as the body is moved through the tube.




Siemens CT Scanner at the University of Virginia

## Prevalence of Hemorrhage in Training Data

7% of training images show hemorrhage of which the most common type is subdural hemorrhage which occurs in 43% of cases






<sup>\*</sup>Images with hemorrhage can have multiple types of hemorrhage

## Predicting Hemorrhage

The model can predict hemorrhage with ~90% accuracy



As of now, our network can correctly classify ~90% of CT scans

Can radiologists take long vacations? Probably not yet, but the tool can be a good double-check or first response when radiologists are not on call

## Uses of the Classifier

The classifier can be used in various situations where cost and speed of hemorrhage detection are a factor









#### **Disaster Response**

Models can often be deployed in a fraction of the time of human resources

#### Rural and High Cost Areas

In areas areas that are currently underserved due to the cost constraints, a model can provide a low cost alternative

#### Radiologist's Assistant

Models can serve as assistants to radiologists by providing a low-cost double check of results

## Catching Fringe Cases

With improved training, the model would have the ability to classify fringe cases that may be hard to predict otherwise

## Future improvements

# Avenues for future improvements include addition of dropout layers as well as use of portions of pre-trained models

#### **Dropout layers to reduce overfitting**

at present we experience overfitting to the training data. Addition of dropout layers

#### **Bagging of various models**

Models of varying performances could be combined to reduce overfitting while increasing accuracy


#### Addition of pre-trained models

Pretrained models such as VGG16 and InceptionV3 have been shown to be effective in CT scan classification

#### **Cost function changes**

In application, different scenarious will have different cost distributions. The model should incorporate this cost function during training

### Sample Training Run Showing Overfitting



## Thank You

Special thanks to the amazing professors and coaches!!



Selecting the Best Density Window


During image preprocessing we will choose the relevant window density from our images that captures blood, bones and brain matter

#### What is a DICOM file?

CT scans used in this model are stored as Digital Imaging and Communications in Medicine (DICOM) files. These files contain an image along with various metadata

## What is windowing?

CT images stored DICOMs allow us to see a wide range of density levels. In our image processing, we will extract only the density window that capture blood, bones and brain



## Model Setup

## The Convolutional Neural Network (CNN) uses various layers to identify key features in the image

**Convolution layers** apply conversions or filters across the image that will help detect features such as horizontal or vertical lines

**Activation layers** allow the model to decide which image features should be considered and which ignored

**Pooling layers** condense the image by extracting the most relevant features of each group of pixels

**Dense layers** condense information into fewer categories that we can use for classification

#### Model Layers

Layer (type)	Output	Shape	Param #
conv2d_18 (Conv2D)	(None,	510, 510, 32)	320
activation_27 (Activation)	(None,	510, 510, 32)	0
max_pooling2d_18 (MaxPooling	(None,	255, 255, 32)	0
conv2d_19 (Conv2D)	(None,	253, 253, 32)	9248
activation_28 (Activation)	(None,	253, 253, 32)	0
max_pooling2d_19 (MaxPooling	(None,	50, 50, 32)	0
flatten_9 (Flatten)	(None,	80000)	0
dense_27 (Dense)	(None,	50)	4000050
activation_29 (Activation)	(None,	50)	0
dense_28 (Dense)	(None,	6)	306
dense_29 (Dense)	(None,	1)	7
Total params: 4,009,931 Trainable params: 4,009,931 Non-trainable params: 0			